


P30114



LOCUS pUG6 4009 bp DNA CIRCULAR SYN 20-JUN-2013
 DEFINITION PCR template vector pUG6, complete sequence.
 ACCESSION P30114
 KEYWORDS .
 SOURCE Unknown.
 ORGANISM Unknown
 Unclassified.
 REFERENCE 1 (bases 1 to 4009)
 AUTHORS Self
 JOURNAL Unpublished.
 COMMENT Gldener, U., Heck, S., Fiedler, T., Beinhauer, J., and Hegemann, J.H. 1996: A new efficient gene disruption cassette for repeated use in budding yeast. Nucleic Acids Research 24, 2519-2524
 FEATURES Location/Qualifiers
 misc_feature 53..86
 /gene="loxP"
 /SECDrawAs="Region"
 misc_feature 96..484
 /gene="pAgTEF"
 /SECDrawAs="Region"
 CDS 485..1294
 /gene="kanMX"

```

        /product="geneticin resistance protein"
        /SECDrawAs="Gene"
misc_feature 1299..1534
        /gene="tAgTEF"
        /SECDrawAs="Region"
misc_feature 1560..1593
        /gene="loxP"
        /SECDrawAs="Region"
misc_feature 1906..2000
        /gene="ori"
        /SECDrawAs="Region"
CDS          complement (2682..3542)
        /gene="AmpR"
        /product="ampicillin resistance protein"
        /SECDrawAs="Gene"
BASE COUNT   1026 a     953 c     949 g     1081 t
ORIGIN
    1 gaacgcgggcc gccagctgaa gtttcgtacg ctgcagggtcg acaaccctta atataacttc
   61 gtataatgta tgctatacga agttattagg tctagagatc tgttttagctt gcctcgtccc
  121 cgccgggtca cccggccagc gacatggagg ccagaatac cctccttgac agtcttgacg
  181 tgcgcagctc aggggcatga tgtgactgtc gcccgtagat ttagcccata catccccatg
  241 tataatcatt tgcattccata ctttttgatg gccgcacggc gcgaagcaaa aattacggct
  301 cctcgttgca gacctgcgag cagggaaaacg ctccccctcac agacgcgttg aattgtcccc
  361 acgccgcgcc cctgtagaga aatataaaaag gttaggattt gccactgagg ttcttctttc
  421 atatacttcc ttttaaaatc ttgctaggat acagttctca catcacatcc gaacataaac
  481 aacctatgggt aaggaaaaga ctacagtttc gaggccgcga ttaaattcca acatggatgc
  541 tgatttatat gggatataaat gggctgcgca taatgtcggg caatcagggtg cgacaatcta
  601 tcgattgtat gggaaagccc atgcgccaga gttgtttctg aaacatggca aaggtagcgt
  661 tgccaatgat gttacagatg agatggtcag actaaactgg ctgacggaat ttatgcctct
  721 tccgaccatc aagcatttta tccgtactcc tgatgatgca tggttactca ccaactgcgat
  781 ccccgcaaaa acagcattcc aggtattaga agaatactct gattcagggtg aaaatattgt
  841 tgatgcgctg gcagtgttcc tgcgcgggtt gcattcgatt cctgtttgta attgtccttt
  901 taacagcgat cgcgtatttc gtctcgctca ggcgcaatca cgaatgaata acggtttggg
  961 tgatgcgagt gattttgatg acgagcgtaa tggctggcct gttgaacaag tctggaaaga
 1021 aatgcataag cttttgccat tctcaccgga ttcagtcgtc actcatgggtg atttctcact
 1081 tgataacctt atttttgacg aggggaaatt aataggttgt attgatgttg gacgagtcgg
 1141 aatcgcagac cgataccagg atcttgccat cctatggaac tgcctcgggtg agttttctcc
 1201 ttcattacag aaacggcttt ttcaaaaata tggattgatg aatcctgata tgaataaatt
 1261 gcagtttcat ttgatgctcg atgagttttt ctaatcagta ctgacaataa aaagattctt
 1321 gttttcaaga acttgtcatt tgtatagttt ttttatattg tagttgttct attttaatca
 1381 aatgttagcg tgatttataat tttttttcgc ctgcacatca tctgccaga tgccaagtta
 1441 agtgcgcaga aagtaatatc atgcgtcaat cgtatgtgaa tgctggtcgc tatactgctg
 1501 tcgattcgat actaacgccg ccattccagtg tcgaaaacga gctctcgaga acccttaata
 1561 taacttcgta taatgtatgc tatacgaagt tattaggtga tatcagatcc actagtggcc
 1621 tatgcggccg cggatctgcc ggtctcccta tagtgagtcg tattaatttc gataagccag
 1681 gtttaacctg attaatgaat cggccaacgc gcggggagag tctggttgcg tattggcgcg
 1741 tcttccgctt cctcgctcac tgactcgctg cgctcggtcg ttcgggtcgc gcgagcggtg
 1801 tcagctcact caaaggcggg aatacggtta tccacagaat caggggataa cgcaggaaaag
 1861 aacatgtgag caaaaggcca gcaaaaggcc aggaaccgta aaaaggccgc gttgctggcg
 1921 tttttccata ggctccgccc ccctgacgag catcacaaaa atcgacgctc aagtcagagg
 1981 tggcgaaacc cgacaggact ataaagatac caggcgtttc cccctggaag ctccctcgtg
 2041 cgctctcctg ttccgacctt gccgcttacc ggatacctgt ccgcctttct cccttcggga
 2101 agcgtggcgc tttctcaatg ctacgctgtg aggtatctca gttcggtgta ggtcgttcgc
 2161 tccaagctgg gctgtgtgca cgaacccccg gttcagcccc accgctgcgc cttatccggt
 2221 aactatcgtc ttgagtccaa cccggtaaga cactgacttat cgccactggc agcagccact
 2281 ggtaacagga ttagcagagc gaggtatgta ggcggtgcta cagagttctt gaagtgggtg
 2341 cctaactacg gctacactag aaggacagta tttggtatct gcgctctgct gaagccagtt
 2401 accttcggaa aaagagtttg tagctcttga tccggcaaac aaaccaccgc tggtagcggg
 2461 ggtttttttg tttgcaagca gcagattacg cgcagaaaaa aaggatctca agaagatcct
 2521 ttgatctttt ctacgggggtc tgacgctcag tggaaacgaaa actcacgtta agggattttg
 2581 gtcattgagat tatcaaaaag gatcttcacc tagatccttt taaattaaaa atgaagtttt
 2641 aaatcaatct aaagtatata tgagtaaact tgggtctgaca gttaccaatg cttaatcagt
 2701 gaggcaccta tctcagcgat ctgtctattt cgttcatcca tagttgcctg actccccgct

```

2761	gtgtagataa	ctacgatacg	ggaggggctta	ccatctggcc	ccagtgtctgc	aatgataccg
2821	cgagacccac	gctcaccggc	tccagattta	tcagcaataa	accagccagc	cggaagggcc
2881	gagcgcagaa	gtggtcctgc	aactttatcc	gcctccatcc	agtctattaa	ttgttgccgg
2941	gaagctagag	taagtagttc	gccagttaat	agtttgcgca	acgttggtgc	cattgctaca
3001	ggcatcgtgg	tgtcacgctc	gtcgtttggg	atggcttcat	tcagctccgg	ttcccaacga
3061	tcaaggcgag	ttacatgata	ccccatggtg	tgcaaaaaag	cggttagctc	cttcgggtcct
3121	ccgatcgttg	tcagaagtaa	gttggccgca	gtgttatcac	tcattggttat	ggcagcactg
3181	cataattctc	ttactgtcat	gccatccgta	agatgctttt	ctgtgactgg	tgagtactca
3241	accaagtcata	tctgagaata	gtgtatgcgg	cgaccgagtt	gctcttgccc	ggcgtcaata
3301	cgggataata	ccgcgccaca	tagcagaact	ttaaaagtgc	tcattcattgg	aaaacgttct
3361	tcggggcgaa	aactctcaag	gatcttaccg	ctgttgagat	ccagttcgat	gtaaccact
3421	cgtgcaccca	actgatcttc	agcatctttt	actttcacca	gcgtttctgg	gtgagcaaaa
3481	acaggaaggc	aaaatgccgc	aaaaaaggga	ataagggcga	cacggaaatg	ttgaatactc
3541	atactcttcc	tttttcaata	ttattgaagc	atttatcagg	gttattgtct	catgagcgga
3601	tacatatattg	aatgtattta	gaaaaataaa	caaatagggg	ttccgcgcac	atttccccga
3661	aaagtgccac	ctgacgtcta	agaaaccatt	attatcatga	cattaacctc	taaaaatagg
3721	cgtatcacga	ggccctttcg	tctcgcgcgt	ttcgggtgatg	acggtgaaaa	cctctgacac
3781	atgcagctcc	cggagacggg	cacagcttgt	ctgtaagcgg	atgccgggag	cagacaagcc
3841	cgtcagggcg	cgtcagcggg	tgttgggcggg	tgtcggggct	ggcttaacta	tgcggcatca
3901	gagcagattg	tactgagagt	gcaccatatg	gacatattgt	cgtagaacg	cggctacaat
3961	taatacataa	ccttatgtat	catacacata	cgatttaggt	gacactata	

//